

# In-depth analysis and optimization of the ProTrap XG sample preparation device using quantitative DIA mass spectrometry.

Jean-Philippe Couture<sup>1</sup>, Jean-François Noël<sup>1</sup>, Victoria Miller<sup>2</sup>, Sara Lahsae Little<sup>2</sup>, and Hugo Gagnon<sup>1</sup>.

1- PhenoSwitch Bioscience, Sherbrooke, QC, Canada. 2- Proteoform Scientific, Halifax, NS, Canada.

## OVERVIEW

The ProTrap XG is an all-in-one, fast, and robust proteomics sample preparation device that performs close to our optimized and fine-tuned protocol.

## MATERIALS AND METHODS

Side-by-side variations in the protocols

PhenoSwitch No Device	Proteoform ProTrap XG	Modified ProTrap XG
No salt during precipitation	50 mM NaCl during precipitation	50 mM NaCl during precipitation
Protein pellet wash with methanol	Protein pellet wash with acetone	Protein pellet wash with methanol
Reduction with 10 mM DTT	Reduction with 20 mM DTT	Reduction with 10 mM DTT
Alkylation with 15 mM IAA	Alkylation with 100 mM IAA	Alkylation with 15 mM IAA
Quench with 10 mM DTT	No quenching	Quench with 10 mM DTT
Digestion in 0.75 M urea + 50 mM Tris pH 8.0	Digestion in 2 mM CaCl <sub>2</sub> + 50 mM Tris pH 8.0	Digestion in 0.8 M urea + 50 mM Tris pH 8.0
30:1 protein:enzyme ratio	100:1 protein:enzyme ratio	30:1 protein:enzyme ratio
Stop digestion with 2% FA	Stop digestion with 0.1% TFA	Stop digestion with 2% FA

Side-by-side variations in the SPE


Steps	ProTrap XG SPE	Phenomenex SPE
Organic priming 1	100% MeOH	100% ACN
Organic priming 2	50% ACN + 0.1% TFA	N/A
Aqueous wash	5% ACN + 0.1% TFA	2% FA
Sample loading	Sample in 5% ACN + 0.1% TFA	Sample in 2% FA
First Wash	5% ACN + 0.1% TFA	0.1% TFA
Second Wash	N/A	2% FA
Elution	50% ACN + 0.1% TFA	40% ACN + 2% FA

### Mass Spectrometry Method

Sample acquisition was performed on a TripleTOF 6600 (Sciex) coupled to an Eksigent MicroLC200. Peptides were separated on a Kinetex XB column (150 x 0.3 mm, 2.6 µm particle) over a 60 minute LC gradient (A: 3% DMSO + 0.2% FA in water, B: 3% DMSO + 0.2% FA in ethanol). Acquisitions were performed in data independent acquisition, with 10 m/z windows ranging from 350 to 1250 m/z. Peptides were quantified using either the SWATH 2.0 micro app (Peakview, Sciex) and DIA-NN (Demichev et. al. Nature methods, 2020).


## RESULTS

**Fig.1:** Normalized peptide signals over a 60 minutes gradient, for the different sample preparation protocols. Each protocol was repeated 4 times.




**Fig.2:** Global comparison of the different sample preparation protocols.

TEST A: No device, Protocol Standard PhenoSwitch, SPE Phenomenex  
TEST B: PTXG Device, PTXG Protocol, PTXG SPE  
TEST C: PTXG Device, PTXG Protocol, Phenomenex SPE  
TEST D: PTXG Device, PTXG Modified Protocol, PTXG SPE  
TEST E: PTXG Device, PTXG Modified Protocol, Phenomenex SPE



**Fig.3:** Heatmap clustering analysis of peptides with different PTMs or miscleavages/missed cleavages, for all different sample preparation protocols.



**Fig.4 :** Head-to-head comparison of the ProTrap XG and PhenoSwitch Bioscience's current optimized protocol.

Comparisons	ProTrap XG	PhenoSwitch
Complete workflow time	2 days	2 days
Hands-on time	2.5 hours	5 hours
Pricing	\$\$	\$\$\$
Results	Comparable	Optimized
Convenience	Everything in one device	Precipitation, digestion, and SPE on 3 different devices

## CONCLUSIONS

- The ProTrap XG device performs similarly to PhenoSwitch Bioscience's current optimized protocol.
- The biggest differences are attributable to the SPE step (more miscleavages/missed cleavages, longer peptides).
- Given the savings in preparation time and resources, we conclude that the ProTrap XG device is a viable sample preparation tool for whole cell proteomics.

